



# GUARDIANPLUS

## Sedation Monitor

The Guardian Plus monitor has an exceptional user-interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Experience the immediate and measurable benefits of the smartphone-like user-interface that delivers 90% of all common monitoring functions with just two steps.



**3.2 KG**

Weight

**4 HR**

Battery  
Life

**8**

Waveforms

**10 IN**

Touch Screen  
Display

# Features

---

- 10" Touchscreen with Gesture Control
- Customizable Alarms and Thresholds
- Adult, Pediatric, and Neonate Presets
- 18 Configurable Short Cuts and Quickkeys
- Integrated Co2 Option

## Standard Parameters:

ECG (3-5 Lead), HR, RESP, Masimo SET SpO2, NIBP, PR, TEMP

## Optional Parameters:

Necollar OxiMax SPO2, CO2

# Clinician-focused Technology for Enhanced Patient Care

---

- Parameter list tile with automatic highlight of abnormal readings
- One touch to drag icon tab right and instantly view minitrend screen or set baseline vitals at the start of a patient case
- Remote view and alarm watch bed-to-bed allows up to 12 patients to be viewed
- Serial DIAP and HL7 communication protocols provide seamless transfer of data to EMR



Primary Care



ASC



Dental





Easy Transport



Pole Mount



Bedside Mount

## Ergonomic Design

---

### Delivers Exceptional Performance and Reliability:

- Manufactured with industry-leading plastics that stand up to the most stringent medical grade cleaners
- IP44 rated, resistant to damage often caused by continuous cleaning and usage
- The concealed grab-and-go handle eliminates hard to clean crevices and eases clinical workflow while on the move
- Flexible mounting for bedside and transport provides options to fit varying care environments

## Unmatched Visual Experience

---

- Smartphone-like Multi-gesture, Capacitive Touchscreen Display
- HD Display With Wide Angle View  $\geq 170^\circ$
- Flat User Interface Design, Common Across Platform
- Reach 90% Of Common Monitor Functions in 2 Steps
- Up to 18 Configurable Quick Keys

# Accessories

---

Accessories include one battery, one adult cuff, 3.5 m NIBP hose, one line cord, one roll of paper, ECG N/T cable, 3-lead adult/pediatric snap lead wire set, Masimo or Nellcor SpO2 cable and sensor, temperature probe adapter cable, quick reference guide, and one operator's manual. **All other accessories must be ordered separately.**

## BATTERIES

- Lithium Ion Battery — *115-065140-00*

## CABLES

- Power Cord, 110 Volt — *0012-25-0001*
- Nurse Call Cable — *009-003116-00*
- Analog Output Cable — *009-003117-00*
- Synchronization Cable — *009-003118-00*
- Serial Port Cable — *009-009117-00*

## ECG TRUNK CABLES

- ECG Cable, 10' (3.1 m), Defibrillation-proof, for N/T — *009-005266-00*
- ECG Cable, 20' (6.2 m), Defibrillation-proof, for N/T — *009-005267-00*
- ECG Cable, 10' (3.1 m), ESU-Proof, for N/T — *009-005268-00*
- ECG Cable, 20' (6.2 m), ESU-Proof, for N/T — *009-005269-00*

## ECG LEADWIRES-3 LEAD PINCH/CLIP

- 3-Lead, N/T, AHA, pinch/clip, 24" — *009-004771-00*
- 3-Lead, N/T, AHA, pinch/clip, 36" — *009-004772-00*

## ECG LEADWIRES-3 LEAD SNAP

- 3-Lead, N/T, AHA, snap, 24" — *009-004765-00*
- 3-Lead, N/T, AHA, snap, 36" — *009-004766-00*

## ECG LEADWIRES-5 LEAD PINCH/CLIP

- 5-Lead, N/T, AHA, pinch/clip, 24" — *009-004786-00*
- 5-Lead, N/T, AHA, pinch/clip, 36" — *009-004787-00*

## ECG LEADWIRES-5 LEAD SNAP

- 5-Lead, N/T, AHA, snap, 24" — *009-004782-00*
- 5-Lead, N/T, AHA, snap, 36" — *009-004783-00*

## PAPER

- Paper, Thermal Strip Local Printer (12/rolls) — *0683-00-0505-02*

# Accessories (Continued)

## TEMPERATURE ACCESSORIES

- Esophageal stethoscope, 400 Series disposable temperature probe, 12 Fr, ES400-12 (20/box) — **0206-03-0112-02**
- Esophageal stethoscope, 400 Series disposable temperature probe, 18 Fr, ES400-18 (20/box) — **0206-03-0118-02**
- Esophageal/rectal, 400 Series disposable temperature probe, 9 Fr, ER400-9 (20/box) — **0206-03-0209-02**
- Esophageal/rectal, 400 Series disposable temperature probe, 12 Fr, ER400-12 (20/box) — **0206-03-0209-02**
- Skin, 400 Series disposable temperature probe, SK 400 (20/box) — **0206-03-0300-02**
- Temperature transition cable for 400 Series disposable probes — **040-001235-00**
- MR401B reusable temperature probe, adult, esophageal/rectal, 2 pin — **0011-30-37392**
- MR403B reusable temperature probe, adult, skin, 2 pin — **0011-30-37393**
- MR402B reusable temperature probe, pediatric/infant, esophageal/rectal, 2 pin — **0011-30-37394**
- MR404B reusable temperature probe, pediatric/infant, skin, 2 pin — **0011-30-37395**
- Temperature probe, disposable, adult/pediatric, esophageal/rectal, 9FR (20/pkg) — **040-003294-00**
- Temperature probe, disposable, all sizes, skin surface (20/pkg) — **040-003295-00**

## SPO2 ACCESSORIES - MASIMO RD SET

### Reusable

- SpO2 sensor, RD SET DCI, Adult (> 30 kg) — **040-003376-00**
- SpO2 sensor, RD SET DCI, Pediatric (10-50 kg) — **040-003377-00**
- SpO2 sensor, RD Set TC-I, tip-clip ear sensor, 3ft, Adult (> 30 kg) — **040-003380-00**

### Disposable

- SpO2 sensor, RD SET Adt, adhesive, Adult (> 30 kg, pkg/20) — **040-003382-00**
- SpO2 sensor, RD SET PDT, adhesive, Pediatric (10 to 50 kg, pkg/20) — **040-003383-00**
- SpO2 sensor, RD Set Inf, adhesive, Infant (3 to 20 kg, pkg/20) — **040-003384-00**
- SpO2 sensor, RD Set Neo, adhesive, Neonatal /Adult (<3kg or >40kg, pkg/20) — **040-003385-00**
- SpO2 sensor, RD Set NeoPt, adhesive, Neonatal (<1 kg pkg/20) — **040-003386-00**
- SpO2 sensor, RD Set NeoPt-500, non-adhesive, (<1 kg pkg/20) — **040-003387-00**

### Cables

- Masimo RD SET cable, 8-pin — **040-003310-00**
- Masimo RD to LNCS adapter cable — **040-003381-00**
- Masimo LNCS to RD adapter — **040-003426-00**

## SPO2 ACCESSORIES, NELLCOR

- Nellcor SpO2 extension cable, 8pin — **0010-20-42712**

## NIBP ACCESSORIES

- Nibp Tubing, Adu/Ped/Neo, with Connectors (3m) — **6200-30-09688**
- Adult Long Disposable Cuff, 25 to 35 cm (Limb) (10/box) — **115-027568-00**
- Large Adult Disposable Cuff, 33 To 47 cm (Limb) (10/box) — **115-027566-00**
- Large Adult Long Disposable Cuff, 33 to 47 cm (Limb) (10/box) — **115-027569-00**
- Adult Thigh Disposable Cuff, 46 to 66 cm (Limb) (5/box) — **115-027567-00**

# Accessories (Continued)

---

## CO2 ACCESSORIES

- Dryline II Water Trap, Adult/Pediatric (10/box) — **115-058733-00**
- Dryline II Water Trap, Neonatal (10/box) — **115-058734-00**
- Sampling Line, Adult/Pediatric, 2.5 M, (25/box) — **115-043017-00**
- Sampling Line, Neonatal, 2.5 M, (25/box) — **115-043018-00**
- Dryline Airway Adapter, Adult/Pediatric, Straight (10/box) — **115-043020-00**
- CO2 Nasal Sample Cannula, Adult, 7' Line (25/box) — **115-043001-00**
- CO2 Nasal Sample Cannula, Pediatric, 7' Line (25/box) — **115-043002-00**
- CO2 Nasal Sample Cannula, Infant, 7' Line (25/box) — **115-043003-00**

## MOUNTING

- Wall Mount with Quick Release - M Series — **045-003427-00**
- M Series 16' Arm (Mounting For Anesthesia) with Quick Release — **121-001937-00**
- VHM Mount with Quick Release — **045-002148-00**
- Rolling Stand with Quick Release Mount, Cable Management, One Storage Basket — **045-002960-00**
- Storage Basket, Additional (2nd) Without Handle — **045-002936-00**
- Cleaning/disinfecting Wipes Mount — **045-002935-00**
- Bedrail Hook — **8000-30-90169**
- Bedrail Clamp — **8000-30-90170**
- Cable Hooks / Attach to M-series Or Vhm Arm — **042-007382-00**
- Quick Release Mount for Arm Mount and Rolling Stand — **045-003424-00**